

MOBILITY



The Lumex Folding Walkers are designed to offer consumers and clinicians a high quality walking aid that gives the user more independence in their everyday lives. The Lumex walkers are constructed using lightweight aluminum making them easy to lift and maneuver. The walkers feature side braces for stability and a secure handhold and are offered in both adult and junior sizes.

Dual-Release Folding Walker

- Raised H-frame design facilitates over-toilet positioning
- Rigid dual side brace design
- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide caps to facilitate use on most surfaces for wheeled options
- Item No. 616070A, 616270A-4, 616370A-4



E0135*

Single-Release Folding Walker

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Raised H-frame design facilitates over-toilet positioning
- Rigid lower side brace design
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Item No. 613070A



E0135*

Everyday Dual-Release Folding Walker

- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide tips to enhance maneuverability on most surfaces for wheeled options
- Item No. 716070A-2, 716070A-4, 716170A-4, 716270A-2, 716270A-4, 716370A-2, 716370A-4



E0135*

Rigid Walker

- Non-folding design provides extra rigidity for use where space is not an issue
- Angled adjustable legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Item No. 1925A



E0135*

Bariatric Dual-Release X-Wide Walker

- Increased width and depth accommodates larger users
- Frame is designed to be robust and rigid to provide additional strength
- Available with and without double 5" fixed wheels
- Item No. 604070A, 604070W (with wheels)



E0135*

colorSelect

ColorSelect Dual-Release Folding Walker

- Same features as the Everyday Dual-Release Walker but available in a variety of color options to appeal to the user
- Item No.* 716270B-2, 716270BK-2, 716270G-2, 716270R-2, 716270P-2, 716270PK-2
- *1 ea also available

