

MOBILITY



The Lumex Folding Walkers are designed to offer consumers and clinicians a high quality walking aid that gives the user more independence in their everyday lives. The Lumex walkers are constructed using lightweight aluminum making them easy to lift and maneuver. The walkers feature side braces for stability and a secure handhold and are offered in both adult and junior sizes.

Dual-Release Folding Walker

- Raised H-frame design facilitates over-toilet positioning
- Rigid dual side brace design
- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide caps to facilitate use on most surfaces for wheeled options
- Item No. 616070A, 616270A-4, 616370A-4



Single-Release Folding Walker

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Raised H-frame design facilitates over-toilet positioning
- Rigid lower side brace design
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Item No. 613070A



Everyday Dual-Release Folding Walker

- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide tips to enhance maneuverability on most surfaces for wheeled options
- Item No. 716070A-2, 716070A-4, 716170A-4, 716270A-2, 716270A-4, 716370A-2, 716370A-4



Rigid Walker

- Non-folding design provides extra rigidity for use where space is not an issue
- Angled adjustable legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Item No. 1925A



Bariatric Dual-Release X-Wide Walker

- Increased width and depth accommodates larger users
- Frame is designed to be robust and rigid to provide additional strength
- Available with and without double 5" fixed wheels
- Item No. 604070A, 604070W (with wheels)



colorSelect

ColorSelect Dual-Release Folding Walker

- Same features as the Everyday Dual-Release Walker but available in a variety of color options to appeal to the user
- Item No.* 716270B-2, 716270BK-2, 716270G-2, 716270R-2, 716270P-2, 716270PK-2
- *1 ea also available



Platform Walker Attachment

- Designed for users with limited strength or dexterity who are unable to grip the walker
- Adjustable strap comfortably secures arm in the forearm rest
- Padded forearm rest can be pivoted and angled to fit user comfortably and to assist in lifting walker
- Universal clamp fits most 1" tubing walkers
- Tool-free installation and adjustment allows quick fitting for most patient sizes and types
- Item No. 6132B



UpRise Onyx Folding Walker

- Combination of folding walker, rising aid and toilet frame all in one
- Secondary handles provide stable assistance from a seated to a standing position
- Attractive, upscale Onyx finish
- Single release folding mechanism
- Durable aluminum tubing provides strength while remaining lightweight
- Adjustable height fits a broad range of user heights
- Item No. 700175C-2, 700175CR (Retail Packaged)



Walkane Hemi Walker

- Lighter than a walker, more stable than a quad cane
- Ideal for users who do not require a walker but need a wider base for support and for individuals with the use of only one hand or arm
- Contoured, foam hand grips are comfortable and provide secure hand hold
- Comfortable, dual level handgrips provide support from a seated to standing position
- Item No. 6015A



*HCPCS Code

Item No.	Width Inside Handgrips	Overall Width	Width Inside Rear Legs	Overall Depth	Minimum Height	Maximum Height	Weight Capacity ¹ EVENLY DISTRIBUTED
616070A	16.5 in / 41.9 cm	22.7 in / 41.9 cm	20.1 in / 51 cm	16.5 in / 41.9 cm	32.5 in / 82.5 cm	41.5 in / 105.4 cm	300 lb / 136 kg
616270A-4	16.5 in / 41.9 cm	26 in / 66 cm	20.1 in / 51 cm	16.5 in / 41.9 cm	34.5 in / 87.6 cm	41.5 in / 105.4 cm	300 lb / 136 kg
616370A-4	16.5 in / 41.9 cm	25.3 in / 64.2 cm	18.6 in / 47.2 cm	15.5 in / 39.3 cm	29.5 in / 74.9 cm	34.5 in / 87.6 cm	300 lb / 136 kg
604070A	20 in / 50.8 cm	26.1 in / 66.2 cm	23 in / 55.4 cm	18 in / 45.7 cm	31 in / 78.7 cm	38 in / 96.5 cm	600 lb / 272.1 kg
604070W	20 in / 50.8 cm	29.2 in / 74.1 cm	23 in / 55.4 cm	19 in / 48.2 cm	33 in / 83.8 cm	40 in / 101.6 cm	600 lb / 272.1 kg
716070A-1 716070A-2 716070A-4	17 in / 43.1 cm	22 in / 55.8 cm	20 in / 50.8 cm	18.5 in / 46.9 cm	32 in / 81.2 cm	38 in / 96.5 cm	300 lb / 136 kg
716170A-4	17 in / 43.1 cm	22 in / 55.8 cm	20 in / 50.8 cm	18.5 in / 46.9 cm	25 in / 63.5 cm	32 in / 81.2 cm	300 lb / 136 kg
716270A-1 716270A-2 716270A-4 716270 ColorSelect series	17 in / 43.1 cm	25.5 in / 64.7 cm	20 in / 50.8 cm	20.2 in / 51.3 cm	32 in / 81.2 cm	39 in / 99 cm	300 lb / 136 kg
716370A-2 716370A-4	17 in / 43.1 cm	24 in / 60.9 cm	20 in / 50.8 cm	18 in / 45.7 cm	25.5 in / 64.7 cm	32.5 in / 82.5 cm	300 lb / 136 kg
613070A	17.25 in / 43.8 cm	23.2 in / 58.9 cm	19.7 in / 50 cm	19.5 in / 49.5 cm	32.5 in / 82.5 cm	39.5 in / 100.3 cm	300 lb / 136 kg
700175C-2 700175CR	18 in / 45.7 cm	26 in / 66 cm	23 in / 55.4 cm	20 in / 50.8 cm	32 in / 81.2 cm	39 in / 99 cm	400 lb / 181.4 kg
1925A	17 in / 43.1 cm	22.5 in / 57.1 cm	19.7 in / 50 cm	20.5 in / 52 cm	31 in / 78.7 cm	38 in / 96.5 cm	400 lb / 181.4 kg
6015A	-	16.5 in / 41.9 cm	-	13.2 in / 33.5 cm	29.5 in / 74.9 cm	34.5 in / 87.6 cm	300 lb / 136 kg