

FOLDING WALKERS

EVERYDAY DUAL RELEASE FOLDING WALKERS

Features

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- Accommodates all leg attachments
- Designed to meet the requirements of HCPCS codes:
(without wheels) E0135
(with wheels) E0143
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty



			5" Fixed Wheels		No wheels		Width inside handgrips		Width inside rear legs		overall Width		Depth front to rear		Depth folded		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
716070A-2	Adult		●				18	45.72	21	53.34	22	55.88	17	43.18	4	10.16	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	2/cs
716070A-4	Adult		●				18	45.72	21	53.34	22	55.88	17	43.18	4	10.16	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs
716170A-4	Junior		●				16.5	41.91	19	48.26	19.5	49.53	15	38.10	4	10.16	25.5	64.77	34.25	87.00	4' 7" - 5' 3"	1.4 - 1.6	300	136.08	4/cs
716270A-2	Adult		●				18	45.72	21	53.34	22	55.88	17	43.18	4	10.16	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	2/cs
716270A-4	Adult		●				18	45.72	21	53.34	22	55.88	17	43.18	4	10.16	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs
716370A-2	Junior		●				16.5	41.91	19	48.26	19.5	49.53	15	38.10	4	10.16	25.5	64.77	34.25	87.00	4' 7" - 5' 3"	1.4 - 1.6	300	136.08	2/cs
716370A-4	Junior		●				16.5	41.91	19	48.26	19.5	49.53	15	38.10	4	10.16	25.5	64.77	34.25	87.00	4' 7" - 5' 3"	1.4 - 1.6	300	136.08	4/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.