

MOBILITY



ROLLATORS

The Lumex rollators are designed to allow for easy maneuverability and comfort to aid in mobility. Each rollator frame is built for durability and features the following:

- Padded seat for added comfort
- Easy to remove padded backrest in standard and contoured style
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Adjustable handle height to accommodate a broad range of user heights
- Comes with wire basket or pouch for storage of personal items



E0143*
E0156

Walkabout Lite Four-Wheel Rollator

- Lightweight aluminum frame weighs only 14 ½ lb
- 6" wheels designed for easy turning and mobility
- Basket for storage can be placed under seat or in front of frame
- Reinforced frame built for durability
- Available in 8 attractive, retail-friendly color options
- Item No. RJ4300B, RJ4300G, RJ4300K, RJ4300L, RJ4300P, RJ4300R, RJ4300AQ, RJ4300CH



E0143*
E0156

Walkabout Wide Four-Wheel Rollator

- Wider, aluminum frame with 18.5 inch width in between handgrips allows for additional patient comfort when positioning and seated for resting
- Wide 15.5 inch (39.3 cm), padded seat provides larger seating surface for comfort
- 6 inch wheels designed for easy turning and mobility
- Item No. RJ4318B, RJ4318R, RJ4318AQ



E0143*
E0156

Walkabout Basic Four-Wheel Rollator

- Lightweight 13 ½" aluminum frame is an affordable solution to mobility yet maintains ease of use for patient
- Convenient black storage pouch for personal items
- 6" wheels designed for easy turning and mobility
- Item No. RJ4900B, RJ4900R

Walkabout Essentials Four-Wheel Rollator

- Durable, cost effective steel frame
- Padded seat for added comfort
- Removable, hinged, padded backrest can be folded for transport
- 6" wheels designed for easy turning and mobility
- Pouch included for storage of personal items
- Item No. RJ5500R, RJ5500K



E0143*
E0156

Walkabout Junior Four-Wheel Rollator

- Smaller frame is ideal for petite or junior users
- 6" wheels designed for easy turning and mobility
- Tool-free removable padded backrest for comfort
- Item No. RJ4301R



E0143*
E0156

*HCPCS Code



colorSelect
Available on select models



E0143*
E0156

Walkabout Deluxe Contour Four-Wheel Rollator

- Contoured, padded backrest offers added comfort
- 8" wheels are ideal for indoor and outdoor use
- Comes with convenient basket placed under seat or in front frame for personal items
- Reinforced frame is built for durability and long term use
- Padded seat gives user comfort when seated
- Item No. RJ4805B, RJ4805K, RJ4805L, RJ4805P, RJ4805R, RJ4805CH

Walkabout Hemi Four-Wheel Rollator

- Lower seat-to-floor height accommodates shorter individuals
- Convenient wire basket for personal items
- 6" wheels designed for easy turning and mobility
- Item No. RJ4302B, RJ4302K, RJ4302R



E0143*
E0156



E0143*
E0156

Set N' Go Height Adjustable Rollator

- Offers the convenience and flexibility of a two-in-one product all rolled into one
- Easy seat height adjustment to be positioned from 18" to 22" in 1" increments
- Comfortable, seamless padded seat with zippered pouch under seat and larger black pouch
- Adjustable handle height with knobs to secure in place
- 6" wheels designed for easy turning and mobility
- Item No. RJ4700B, RJ4700R, RJ4700S, RJ4700P

Set N' Go Wide Height Adjustable Rollator

- Offers the same features as the standard Set N' Go but in a Wide option with 19" width in between handgrips
- Easy height adjustment on the frame allows seat height to be positioned from 19.5" to 23.5" in 1" increments
- Wider 15.75" comfortable, seamless padded seat with zippered pouch under seat for security of personal items
- Also features Velcro® folding strap to keep rollator folded for storage and transport
- Item No. RJ4718B, RJ4718R, RJ4718S



E0143*
E0156

Walkabout Imperial Bariatric Four-Wheel Rollator

- Black 8" wheels are ideal for indoor and outdoor use
- Padded oversized 20" seat
- Wide frame made of high quality aluminum to remain lightweight
- Comes with large pouch under the seat for storage of personal items
- Easy to remove padded backrest for comfort
- Item No. RJ4400R, RJ4400K (Standard Backbar)
- Item No. RJ4405B, RJ4405K, RJ4405R (Contoured Backbar)

Walkabout Imperial Hemi Bariatric Four-Wheel Rollator

- Same features as the standard Bariatric but in a hemi 19" seat-to-floor height for shorter individuals
- Item No. RJ4402R



E0149*
E0156



Three-Wheel Cruiser

- Provides convenient mobility for those not needing the extra stability of a 4 wheel rollator
- Lightweight aluminum frame
- Available in 3 attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- Item No. 609201B, 609201M, 609201P

E0143*

HybridLX Rollator/Transport Chair

- Lightweight, durable aluminum frame weighs only 22 lbs (this includes footrests)
- Hybrid concept – designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- Ultra-Compact design – easy, tool-free assembly and folds effortlessly for convenient transport and storage
- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling on smooth and other terrain
- Item No. LX1000T, LX1000B



E0143*
E0156



SureStride Rollator

- Patented seat slides along frame allowing user to walk inside the rollator frame allowing for improved gait
- Adjustable push bar handle provides a place to propel rollator forward with better posture and stability
- Ergonomic handgrips with easy brake access
- Intuitive one brake handle can activate both brakes
- Deluxe, soft pouch for personal items can be easily installed and removed
- Large wheels allow for easy maneuverability
- Item No. S9000TS

Item No.	Width Inside Handgrips	Overall Width	Overall Depth	Seat Depth	Seat Width	Seat Height	Handle Height Min/Max	Weight Capacity ¹ EVENLY DISTRIBUTED
RJ4300B, G,K,L,P,R,AQ,CH	16 in / 40.6 cm	24.5 in / 60.9 cm	27.5 in / 69.8 cm	14.5 in / 36.8 cm	13.5 in / 34.2 cm	23 in / 55.4 cm	32 in - 37 in 81.2 cm - 93.9 cm	300 lb / 136 kg
RJ4318B, R,AQ	18.5 in / 46.9 cm	27.5 in / 69.8 cm	27.5 in / 69.8 cm	15.5 in / 39.3 cm	13 in / 33 cm	22.5 in / 57.1 cm	32 in - 37 in 81.2 cm - 93.9 cm	350 lb / 158.7 kg
RJ4900B,R	16.5 in / 41.9 cm	24.2 in / 61.4 cm	24.5 in / 62.2 cm	13 in / 33 cm	12.5 in / 31.7 cm	21 in / 53.3 cm	32.5 in - 37.5 in 82.5 cm - 95.2 cm	250 lb / 113.3 kg
RJ4805B, K,L,P,R,CH	16 in / 40.6 cm	25 in / 63.5 cm	27.5 in / 69.8 cm	14 in / 35.5 cm	13.5 in / 34.2 cm	22.5 in / 57.1 cm	34 in - 38.5 in 86.3 cm - 97.7 cm	300 lb / 136 kg
RJ4301R	12 in / 33.4 cm	20 in / 50.8 cm	24 in / 60.9 cm	12 in / 33.4 cm	12 in / 33.4 cm	18.5 in / 46.9 cm	27.5 in - 32 in 69.8 cm - 81.2 cm	300 lb / 136 kg
RJ4302B,K,R	15.5 in / 39.3 cm	22.5 in / 57.1 cm	21.5 in / 54.6 cm	14 in / 35.5 cm	13.5 in / 34.2 cm	19.5 in / 49.5 cm	32 in - 36.5 in 81.2 cm - 92.7 cm	300 lb / 136 kg
RJ4700B,R,S	16.3 in / 41.4 cm	24.5 in / 60.9 cm	29 in / 73.6 cm	14 in / 35.5 cm	12 in / 33.4 cm	18" - 22"	29 in - 38 in 73.6 cm - 96.5 cm	300 lb / 136 kg
RJ4718B,R,S	19.4 in / 49.2 cm	26.9 in / 68.3 cm	27.5 in / 31.5 in 69.8 cm / 80 cm	12 in / 33.4 cm	15.7 in / 39.8 cm	19.5" - 23.5"	32 in - 38 in 81.2 cm - 96.5 cm	350 lb / 158.7 kg
RJ4400K,R, RJ4405B,K,R	24 in / 60.9 cm	29.5 in / 74.9 cm	27.5 in / 69.8 cm	20 in / 50.8 cm	13.5 in / 34.2 cm	23.5 in / 59.6 cm	34 in - 39 in 86.3 cm - 99 cm	450 lb / 204.1 kg
RJ4402R	22.5 in / 57.1 cm	29.5 in / 74.9 cm	30 in / 76.2 cm	20 in / 50.8 cm	13.5 in / 34.2 cm	20.5 in / 52 cm	31 in - 34 in 78.7 cm - 86.3 cm	450 lb / 204.1 kg
LX1000B,T	18.5 in / 46.9 cm	26 in / 66 cm	40 in / 101.6 cm	28.5 in / 52.1 cm	18 in / 45.7 cm	300 lb / 136 kg	30 in - 37 in 76.2 cm - 93.9 cm	300 lb / 136 kg
S9000TS	16 in / 40.6 cm	24.5 in / 60.9 cm	31.7 in / 80.5 cm	15.1 in / 38.4 cm	17.2 in / 43.7 cm	22 in / 55.9 cm	29.4 in - 35.4 in 74.7 cm - 89.9 cm	300 lb / 136 kg
609201B,M,P	14.5 in / 36.8 cm	25 in / 63.5 cm	22.5 in / 57.1 cm	-	-	-	33 in - 37 in 83.8 cm - 93.9 cm	250 lb / 113.3 kg

Color Select is a registered trademark of the Rollex Corporation

¹ maximum weight capacity - evenly distributed - refer to user manual *HCPCS Code

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

published **June 2026**

Images are digital reproductions. perceived texture and coloring may vary on
finished product. Monitor and printer settings will effect appearance of color.