







Item  
POWDER-COATED METAL  
SAS-TEE4 TEE PROBE STORAGE CABINET  
Locking Overlay Hinged Door with Glass Insert & Sloped Top



120V Hospital Grade  
Plug Exits top right  
hand side with  
approx. 8' of cord

MU-TRUHEPA  
Forced HEPA  
Filtered Air

8x  
MU-RACK-TEE2-SPLIT  
Independently  
Adjustable TEE Probe  
and Transducer Holders  
with Cord Management



003-00-0044  
Padded Back Wall  
Removable  
Drip Tray  
Leveling Feet  
adjustment of ~1/2"



Item  
POWDER COATED METAL  
SAS-TEE8 TEE PROBE STORAGE CABINET  
Locking Overlay Hinged Doors with Glass Inserts & Sloped Top