


Digi-Flex® Hand Exercisers

Select the hand exerciser that's best for you!

- Use to rehabilitate injuries or conditions affecting the hand and upper extremities
- Develops isolated finger strength, flexibility, and coordination
- Available in eight (8) strength levels
- Comes with exercise pamphlet

Digi-Flex®

- **Digi-Flex®** model has upper and lower spring to allow for more hand motion flexibility. Use for individual finger or hand



Digi-Flex LiTE®

- **Digi-Flex LiTE®** is An inexpensive option for patient giveaway or sale. It has individual finger locks and a single spring design



Digi-Flex Thumb®

- **Digi-Flex Thumb®** allows for thumb involvement. The thumb piece can attach at three different positions to accommodate all thumb sizes and swivels for all thumb positions



force to compress
finger/hand (lb)


tan		xx-light	0.75	2.5			
yellow		x-light	1.5	5.0			
red		light	3.0	10.0			
green		medium	5.0	16.0			
blue		heavy	7.0	23.0			
black		x-heavy	9.0	30.0			
silver		xx-heavy	11.0	36.0			
gold		xxx-heavy	13.0	43.0			
yel	red	grn	blu	blk	set of 5 (1 each)		
yel	red	grn	blu	blk	set with plastic stand		
yel	red	grn	blu	blk	set with metal stand		
tan	yel	red	grn	blu	blk	sil	gld
set of 8 (1 each)							
tan	yel	red	grn	blu	blk	sil	gld
set of 8 with 2 metal stands							

Digi-Flex®

Digi-Flex LiTE®

Digi-Flex Thumb®

10-3770	10-3760
10-3771	10-3761
10-3772	10-3762
10-3773	10-3763
10-3774	10-3764
10-3775	10-3765
10-3776	10-3766
10-3777	10-3767
10-3797	10-3769
10-3798	10-3785
10-3799	10-3786
10-3778	10-3768
10-3779	10-3787



Digi-Flex® Stands

- Holds up to 5 Digi-Flex® units

metal stand	10-0750
plastic stand	10-0747

Digi-Flex Multi® Hand Exerciser

Create a customized hand exerciser in 3 easy steps!



pinch release tabs and
pull up to remove button



insert button until
a "snap" is heard



use your customized
hand exerciser

Clinic Sets

Customize hand therapy like never before. Adjust each finger's resistance to fit the patient's rehab program.

Clinic sets are available with either 5 bases or 5 pre-built exercisers, 20 or 32 piece resistive button sets with storage case, and rack.

More bases, exercisers, buttons, and button sets are available for purchase as needed. See the following link for more information: Fab-Ent.com/DigiFlexMulti

small clinic set with 5 bases

10-3840 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3841 standard, 20 buttons 4 ea yel red grn blu blk

large clinic set with 5 pre-built exercisers

10-3845 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3846 standard, 20 buttons 4 ea yel red grn blu blk



small clinic set
with 5 bases,
deluxe (10-3840)



patented



large clinic set with
5 pre-built exercisers,
standard (10-3846)



Individual Finger Locks

Insert the white finger lock button when finger movement or exertion is counter-indicated.



10-3738-4 4 each with case