


NEW

Digi-Flex® LiTE® exerciser

inexpensive enough for patient give-aways
finger locks prevent compensation and allow for finger isolation




- individual finger locks
- single spring design


		force to compress finger / hand (lb)		
tan	xx-light	.75 / 3.0		10-3770
yellow	x-light	1.5 / 6.0		10-3771
red	light	3.0 / 12.0		10-3772
green	medium	5.0 / 20.0		10-3773
blue	heavy	7.0 / 28.0		10-3774
black	x-heavy	9.0 / 36.0		10-3775
silver	xx-heavy	11.0 / 44.0		10-3776
gold	xxx-heavy	13.0 / 52.0		10-3777
yel red grn blu blk	set (1 ea)			10-3797
yel red grn blu blk	set with plastic stand			10-3798
yel red grn blu blk	set with metal stand			10-3799
tan yel red grn blu blk sil gld	set of 8			10-3778
tan yel red grn blu blk sil gld	set with 2 metal stands			10-3779

NEW

Digi-Flex® Thumb™ exerciser



- exercise fingers and thumb
- design allows for both left and right thumb usage
- thumb can attach in 3 positions
- thumb rotates in all planes
- individually locking fingers and thumb
- single spring design



		force to compress finger / hand (lb)		
tan	xx-light	.75 / 3.0		10-3760
yellow	x-light	1.5 / 6.0		10-3761
red	light	3.0 / 12.0		10-3762
green	medium	5.0 / 20.0		10-3763
blue	heavy	7.0 / 28.0		10-3764
black	x-heavy	9.0 / 36.0		10-3765
silver	xx-heavy	11.0 / 44.0		10-3766
gold	xxx-heavy	13.0 / 52.0		10-3767
yel red grn blu blk	set (1 ea)			10-3769
yel red grn blu blk	set with plastic stand			10-3785
yel red grn blu blk	set with metal stand			10-3786
tan yel red grn blu blk sil gld	set of 8			10-3768
tan yel red grn blu blk sil gld	set with 2 metal stands			10-3787

display / storage rack

- holds up to 5 units of:
Digi-Flex®, Digi-Flex® LiTE®,
Digi-Flex® Multi®, Via™,
VariGrip®, Pro™, Digi-Extend®
and Digi-Flex® Thumb™

plastic stand only	10-0747
metal stand only	10-0750

