

PRODUCT DATASHEET

Model: **ER** (see below for details)

Description: **Easy Reader Radiopaque Ruler**

Revision: August 2019

PRODUCT ATTRIBUTES:

- Hanging application
- Lead-free, environmentally friendly, safe for patient contact
- Clearly visible, high-contrast, easy to read radiopaque reference marks
- Made from resilient material that is both flexible and strong
- Available in various lengths with and without mounting bracket

MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	Lexan	98
Filler	Epoxy	1
Text	Tungsten	1

DIMENSIONS AND PERFORMANCE

Product Dimension: See below

Part	Description	Dimensions	Box Dim's & Weight (lb)
ER66SPH	SPH	1747 mm x 40 mm x 100 mm	4x4x73" - 3.75 lb
ER34B	POSPIT2	1010 mm x 40 mm x 100 mm	4x4x44" - 3.50 lb
ER34BT	2 Panel Trans-Former	1010 mm x 40 mm x 100 mm	5x5x44" - 3.50 lb
ER51B	C-GL, POSPGL, Cuattro One	1365 mm x 40 mm x 100 mm	4x4x58" - 3.75 lb
ER51BT	3 Panel Trans-Former	1365 mm x 40 mm x 100 mm	5x5x58" - 3.75 lb
ER34	Any - Tape On	960 mm x 40 mm x 3.2 mm	4x4x42" - 2.50 lb
ER51	Any - Tape On	1315 mm x 40 mm x 3.2 mm	4x4x56" - 2.75 lb

FLAMMABILITY RATING

The product is non-flammable.
UL Rating: UL94V-2.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:
Reina Imaging

DEVICE FAMILY CLASSIFICATION

Country	Class
USA	Class 1 per MDD 93/42/EEC, Annex V, NIST Compliant
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

ReinaImaging