

WEIGHT DISTRIBUTION CAPS

Product Advantages

Lightweight (6lb.)

Small footprint


Slip-resistant

Cutout for tether




Durable

Rated up to 750lb.

Available with clear or opaque top



Standing feet weight bearing exams require the patient to stand on the imaging receptor. Protection to expensive DR panels is paramount. Reina Imaging's new [14x17 Weight Distribution Cap](#) is now lighter, less money at yet maintains the same performance. It is rated up to 750 lbs and now is available with a clear top too.



Technical Data

Outside Dimensions

50.8 x 43.18 x 5.7 cm

Weight Capacity

340 kg