

ROLLATORS

3-WHEEL CRUISER

The GF 3-Wheel Cruiser is designed to provide convenient mobility for those individuals that do not require the extra stability of a 4-wheel rollator. This new light weight aluminum walking aid allows the user to have more freedom in movement to perform daily activities.

Features

- Light weight aluminum frame
- Available in 3 new attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Ergonomic hand grips with locking loop brakes
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- Instructions included
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity
250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Width inside handgrips		Width inside rear legs		overall Width		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
609201B	Bondi Blue	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
609201M	Metallic Burgundy	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
609201P	Pewter Metallic	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
	Carry Pouch					14.17	35.99							22	9.98	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.