

ROLLATORS

WALKABOUT HEMI FOUR-WHEEL ROLLATOR

The Lumex Walkabout Four-Wheel Hemi Rollator is similar to the Walkabout Lite as it also offers a comfortable, lightweight and affordable solution to mobility. The Walkabout Hemi is ideal for the more petite or shorter individual. Designed with an aluminum frame, this lightweight rollator weighs 14 lbs and supports up to a 300 lb weight capacity - **EVENLY DISTRIBUTED**. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Features

- Designed with a lower seat to floor height for shorter users
- Ergonomic hand grips
- Easy to operate locking and quick stop loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Rear reflectors to increase visibility at night
- Removable padded straight backbar
- Adjustable handle height to accommodate broad range of user heights
- 6" wheel for indoor and outdoor use
- Basket included that can be placed underneath seat or in front for carrying of personal items
- Includes easy to read, retail-friendly user manual
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4302B	Royal Blue	17	43.18	23.5	59.69	21.5	54.61	14	35.56	13.5	34.29	20	50.80
RJ4302K	Black	17	43.18	23.5	59.69	21.5	54.61	14	35.56	13.5	34.29	20	50.80
RJ4302R	Burgundy	17	43.18	23.5	59.69	21.5	54.61	14	35.56	13.5	34.29	20	50.80

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4302B	Royal Blue	6	15.24	32	81.28	36.5	92.71	14	6.35	300	136.08	1 ea
RJ4302K	Black	6	15.24	32	81.28	36.5	92.71	14	6.35	300	136.08	1 ea
RJ4302R	Burgundy	6	15.24	32	81.28	36.5	92.71	14	6.35	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.