


ROLLATORS

SET N' GO HEIGHT ADJUSTABLE ROLLATOR

The Lumex Set n' Go Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one. Designed with multiple height settings on the frame, the Set n' Go allows the seat height to be adjusted from a hemi 18" to a standard 22" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

Features

- Easy height adjustment on the frame allows seat height to be adjusted from 18" to 22" in 1" increments
- Comfortable, wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Ergonomic handgrips offer a comfortable grip
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Removable padded backbar with easy release pull knobs
- 6" wheel for both indoor and outdoor use
- Designed to meet the requirements of HCPCS: E0143 and E0156
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



		overall Width		Depth at minimum seat Height		Depth at maximum seat Height		Width inside handgrips		minimum Handle Height		maximum Handle Height		Wheel diameter	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4700B	Blue	24.5	62.23	25	63.50	29	73.66	16.3	41.40	29	73.66	38	96.52	6	15.24
RJ4700R	Burgundy	24.5	62.23	25	63.50	29	73.66	16.3	41.40	29	73.66	38	96.52	6	15.24
RJ4700S	Silver	24.5	62.23	25	63.50	29	73.66	16.3	41.40	29	73.66	38	96.52	6	15.24
RJ4700P	Pink	24.5	62.23	25	63.50	29	73.66	16.3	41.40	29	73.66	38	96.52	6	15.24

		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4700B	Blue	14	35.56	12	30.48	18	45.72	22	55.88	14	6.35	300	136.08	1 ea
RJ4700R	Burgundy	14	35.56	12	30.48	18	45.72	22	55.88	14	6.35	300	136.08	1 ea
RJ4700S	Silver	14	35.56	12	30.48	18	45.72	22	55.88	14	6.35	300	136.08	1 ea
RJ4700P	Pink	14	35.56	12	30.48	18	45.72	22	55.88	14	6.35	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.