



## ROLLATORS

# WALKABOUT BASIC FOUR-WHEEL ROLLATOR

The Lumex Walkabout Basic Four-Wheel Rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 13.5 lbs and supports up to a 250 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

### Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Maximum Weight Capacity  
250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



|         |      | Width<br>inside handgrips |       | overall<br>Width |       | overall<br>Depth |       | Seat Width |       | Seat Depth |       | Seat Height |       |
|---------|------|---------------------------|-------|------------------|-------|------------------|-------|------------|-------|------------|-------|-------------|-------|
|         |      | inch                      | cm    | inch             | cm    | inch             | cm    | inch       | cm    | inch       | cm    | inch        | cm    |
| RJ4900B | Blue | 16.5                      | 41.91 | 24.2             | 61.47 | 24.5             | 62.23 | 13         | 33.02 | 12.5       | 31.75 | 21          | 53.34 |
| RJ4900R | Red  | 16.5                      | 41.91 | 24.2             | 61.47 | 24.5             | 62.23 | 13         | 33.02 | 12.5       | 31.75 | 21          | 53.34 |

|         |      | Wheel<br>diameter |       | minimum<br>Handle Height |       | maximum<br>Handle Height |       | Product<br>Weight |      | Weight<br>Capacity<br>EVENLY DISTRIBUTED |        | UOM  |
|---------|------|-------------------|-------|--------------------------|-------|--------------------------|-------|-------------------|------|------------------------------------------|--------|------|
|         |      | inch              | cm    | inch                     | cm    | inch                     | cm    | lb                | Kg   | lb                                       | Kg     |      |
| RJ4900B | Blue | 6                 | 15.24 | 32.5                     | 82.55 | 37.5                     | 95.25 | 13.5              | 6.12 | 250                                      | 113.40 | 1 ea |
| RJ4900R | Red  | 6                 | 15.24 | 32.5                     | 82.55 | 37.5                     | 95.25 | 13.5              | 6.12 | 250                                      | 113.40 | 1 ea |



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Rollators are not intended to be propelled while seated.  
Do NOT use this rollator as a wheelchair or transport device.