

574G

LUMEX®
by graham-field

Three Position Preferred Care Recliner is the recliner that started it all and a standard in the industry! This first recliner is specifically designed to improve the comfort and mobility of residents in extended care facilities. Ergonomically designed seat and back with varied density foam construction provides correct support. Articulating headrest quickly adjusts to maximize individual comfort requirements.



EXTENDED
CARE



HOME
CARE



**250
POUND™**
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
☒	☐	☐	☐	☐	☐	☒	☒	☒	☒	☒	☐	☐	☐
AREAS OF USE													

HEALTHCARE SEATING SOLUTIONS

graham
field™

Model

574Gxxx

‡ xxx REPRESENTS FABRIC COLOR

Positions: THREE POSITION

Three Positions for Patient Comfort: Upright, First Recline and Second Recline
Recline Position Lock: Allows Caregiver to Secure the Chair in Any of the Three Positions



Features:

- Ergonomically designed seat and back with varied density foam construction provides correct support
- Flex-o-Lator seat support system increases comfort by reducing the risk of "bottoming out"
- Articulating headrest quickly adjusts to maximize individual comfort requirements
- 5 inch casters for easy transport: 2 Rear Swivel and Wheel Locking; 2 Fixed Front Wheels
- Includes activity table that can be fixed in position or left unsecured for easy exit from the recliner
- Durable epoxy coated steel frame is easy to clean and maintain



INCLUDES ACTIVITY TABLE



ARTICULATING HEADREST



5 INCH CASTERS



REAR PUSH HANDLE



ACTIVITY TRAY STORAGE

Accessories:

★ sold separately

item no.†

Headrest Pillow	FR5659xxxJS
Headrest Skirt	HRSxxx
Activity Tray Table	included
IV Pole	2101A
IV Mount	5782RG

Specifications:

Overall Width		Overall Depth		Overall Depth <i>fully reclined</i>		Overall Height		Arm Width		Width <i>between arms</i>		Arm Height <i>from floor</i>		Seat Width		Seat Depth		Seat Height		Back Width		Back Height <i>above seat</i>		Caster Diameter		upright legrest fully rec Weight ¹ Capacity EVENLY DISTRIBUTED				
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	Back Angle	lb	kg†		
24.50	62.23	41.00	104.14	65.00	165.10	47.00	119.38	2.00	5.08	20.00	50.80	26.00	66.04	19.00	48.26	19.00	48.26	20.00	50.80	19.50	49.53	29.00	73.66	5.00	12.70	105	105	110	250	113.40

*Ships Freight: request additional shipping quote

‡ xxx in the product number is the placeholder for fabric choice

† standard rounding method used


¹ maximum weight capacity - evenly distributed - refer to user manual

★ accessories sold separately

Available Colors:



427 BLUE RIDGE



851 DOE SKIN



857 JADE



863 ROSEWOOD




401 VINTAGE GOLD



409 WARM TAUPE



454 ROYAL BLUE



432 IMPERIAL BLUE



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham
field

GF0900049revF26