

AXIGUARD

PRODUCT INFORMATION

COPYRIGHT ©2026. ALL RIGHTS RESERVED



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION

AXIGUARD

RA-AXIGD-PLF50



The AxiGuard provides targeted radiation protection for the underarm area, helping shield a critical region that may be left exposed by standard lead aprons. Worn comfortably beneath an existing radiation apron, AxiGuard adds an extra layer of protection without restricting movement, making it ideal for healthcare professionals performing radiation-intensive procedures.

Designed to fit comfortably under most lead aprons, AxiGuard features a front Velcro closure that allows for quick, easy adjustment and a secure, personalized fit. Available in seven sizes ranging from XS to 3XL, its lightweight construction provides dependable axillary coverage while maintaining freedom of movement throughout the workday.

Suitable for both men and women, the AxiGuard is available in Axion Core™ 3000. It offers 0.50 mm Pb lead equivalency for reliable radiation protection in high-exposure environments.

SPECIFICATIONS

SHIPPING WEIGHT: 144 oz	FEATURES: Lightweight Design, Adjustable Fit, Targeted Axillary Protection, Comfortable Under-Apron Wear, Maintains Freedom of Movement
DIMENSIONS: 24 × 4 × 4 in	CLOSURE: Front Velcro Closure
CORE MATERIAL: Axion Core™ 3000	FIT: Unisex
PROTECTION: 0.50 mm Pb Lead Equivalency	

AXIGUARD SIZE CHART	
Measure around the widest part of your chest	
SIZE	INCHES
XS	29" – 32"
S	33" – 36"
M	37" – 39"
L	40" – 43"
XL	44" – 48"
2XL	49" – 51"
3XL	52" – 55"