

MODEL DH25-49SD

LAB AND PHARMACY FREEZER -15°C TO -25°C

Capacity [Cu. Ft.]	46 Cu. Ft.
Capacity [Liters]	1387 Liters
Temperature Range [°C]	-15°C to -25°C
Temperature Range [°F]	5°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	52" x 33 1/2" x 81"
Interior [Inches]	49" x 26" x 55"
Exterior [Centimeters]	132 x 81 x 206
Interior [Centimeters]	120 x 72 x 140

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Adaptive defrost
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Doors	Two Doors, Keyed Lock
Door Material	Insulated Steel
Interior Light	LED
Porthole	1/2" Diameter Left Side
Casters	4, Locking
Storage	10 Adjustable Shelves
Shipping Weight	610 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	20 Amp Dedicated 13.9A
NEMA Configuration	NEMA 5-20
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Alarm Relay	Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

