

MODEL DHH20-14SDFMS

FLAMMABLE MATERIAL STORAGE FREEZER -15°C TO -25°C

Capacity [Cu. Ft.]	14 Cu. Ft.
Capacity [Liters]	396 Liters
Temperature Range [°C]	-15°C to -25°C
Temperature Range [°F]	5°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	23.75" x 24.75" x 73"
Interior [Inches]	19" x 17.5" x 63.75"
Exterior [Centimeters]	60 x 63 x 185
Interior [Centimeters]	48 x 44 x 162

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Porthole	3/4" Diameter
Leveling Legs	4
Storage	6 Fixed Shelves
Shipping Weight	195 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3.5A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.

Warranty 2 Years Parts & Labor.