

advantage™ plus 6121

Automatic Digital Home Blood Pressure Monitor

ENHANCED VISUALIZATION FOR GREATER INSIGHTS

The ADC Advantage™ Plus 6121 takes home blood pressure monitoring to the next level, with advanced features and deeper insights. Designed for those who want more control and understanding of their readings, it offers ADC's innovative MAM Advantage™ triple-measurement averaging, a color-coded hypertension indicator, and a backlit display for easy visibility in any light. Choose one of three models for a custom cuff fit.




Three Cuff Size Options:

Model 6121: Adult, Extended Range 8.7 to 16.5 in. (22 to 42 cm)
Model 6121SA: Small Adult 6.7 to 8.7 in. (17 to 22 cm)
Model 6121X: Large Adult 12.5 to 20.5 in. (32 to 52 cm)



5-Year
Warranty



Premium Comfort Contour™ Cuffs

At the heart of our upper-arm monitors is our thoughtfully engineered Comfort Contour™ cuff, crafted from soft, high-quality nylon and ergonomically contoured to fit the natural shape of your arm. Unlike flat, straight cuffs, a curved or conical cuff helps distribute pressure more evenly along the arm's contour, reducing pinching or bunching and ensuring the bladder inflates consistently over the arm for a better fit and more accurate readings.

- ✓ **Accurate, Comfortable Readings**
Smart Logic Modulation (SLM) Technology adjusts cuff inflation for a smoother, gentler experience while optimizing for comfort and accuracy.
- ✓ **Automatic Fit Check**
Built-in fit detection ensures your cuff is positioned correctly before each reading, reducing user error.
- ✓ **MAM Advantage™ Averaging Mode**
A technique recommended by the American Heart Association® for enhanced accuracy, our automated averaging mode takes three consecutive measurements then calculates and displays a single, result.
- ✓ **Two-User Memory**
Stores up to 99 readings per user with date and time stamps, making it easy to track trends over time.
- ✓ **Traffic Light Indicator**
This built-in color-coded scale visually aligns readings with the American Heart Association® guidelines for quick interpretation.
- ✓ **Irregular Heartbeat Detection**
Alerts you to potential irregularities so you can share meaningful information with your healthcare provider.
- ✓ **Easy-to-Read Display**
Large, high-contrast screen with oversized numerals simplifies viewing for all users.