

MODEL DHH20-17SDFX

EXPLOSION PROOF FREEZER -15°C TO -25°C

Capacity [Cu. Ft.]	17 Cu. Ft.
Capacity [Liters]	481 Liters
Temperature Range [°C]	-15°C to -25°C
Temperature Range [°F]	5°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	30.75" x 28.5" x 67.75"
Interior [Inches]	26.25" x 22.5" x 58.25"
Exterior [Centimeters]	78 x 72 x 172
Interior [Centimeters]	67 x 57 x 148

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Leveling Legs	4
Storage	6 Fixed Shelves
Shipping Weight	220 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
Electrical Plug	Not Supplied (Hardwired)
Certification(s)	ETL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations
Class 1, Division II, Group C & D**

Warranty 2 Years Parts & Labor.