

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY


SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



LABORATORY RESEARCH & TESTING | PHARMACY STORAGE | GENERAL PURPOSE REFRIGERATION

Laboratory Sample Storage

Our complete line of laboratory and pharmacy freezers provide the reliability and security essential for the protection of sensitive lab samples in the -10°C to -25°C range. Ideal for the demanding requirements of daily lab applications, our range of freezers offers the dependability and safety required for your daily lab applications.

General Purpose Use

So-Low's Lab and Pharmacy freezers were designed specifically for use in a laboratory and pharmacy setting. These freezers offer many storage capabilities such as shelves and racks to store samples, chemicals, and other laboratory products.

Additional Security Options

So-Low Lab and Pharmacy freezers also come with various options, such as dataloggers, chart recorders and audible visual alarm systems for temperature security.

Secure Protection

Units come standard with a door lock, protecting your valuable samples from tampering.

MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
-----------------	---------------	--------	----------------------	------------------------------------	------------------------------------	-----------------	----------------

UNDERCOUNTER

MV20-1.5UCF	1.5	42	-20°C to -20°C	14" x 13.25" x 13.75"	18.5" x 19.5" x 19.5"	MANUAL	30 LBS
MV30-4UCF	3.2	90	-10°C to -25°C	19.5" x 19.25" x 20.75"	27.75" x 27.75" x 34"	MANUAL	130 LBS
MV40-4UCF	3.2	90	-20°C to -40°C	19.5" x 19.25" x 20.75"	27.75" x 27.75" x 34"	MANUAL	130 LBS

GENERAL PURPOSE UPRIGHT

DHK25-10SD	9.53	270	-10°C to -25°C	19.5" x 18" x 48.5"	27.75" x 25.25" x 70.75"	MANUAL	225 LBS
DHW20-30	30	849	-20°C to -20°C	32" x 22.5" x 65"	35.5" x 35" x 72"	MANUAL	350 LBS

SCIENTIFIC UPRIGHT

DHK20-20MDP	20	566	-10°C to -25°C	25.5" x 21" x 64.75"	30.5" x 29.5" x 74.75"	MANUAL	225 LBS
DHW20-30MDP	30	849	-10°C to -20°C	32" x 22.5" x 65"	35.75" x 34.25" x 72.5"	MANUAL	450 LBS



MICROBIOLOGICAL RESEARCH & TESTING | DIVERSE LAB APPLICATIONS | CULTURE STORAGE

High Performance

General purpose incubators deliver a degree of temperature uniformity usually only found in more expensive, application specific models. This high level of performance is due to a precise microprocessor controller and a unique warm air jacket design.

CO2 Incubators

The CO2 Incubator family incorporates more available options than ever before. These CO2 Incubators combine contamination control, decontamination, and precise temperature and CO2 control for years of constant performance.

MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	UNIT WEIGHT
-----------------	---------------	--------	----------------------	------------------------------------	------------------------------------	----------------

GENERAL PURPOSE INCUBATOR

SMI2	2	56	8°C above ambient to 70°C	15" x 15" x 15.25"	21.5" x 22.5" x 26"	100 LBS
SMI6	5.9	166	8°C above ambient to 70°C	19.5" x 20" x 26"	25.5" x 27.5" x 38"	125 LBS
SMI7	6.5	184	8°C above ambient to 70°C	23.75" x 24" x 20"	30" x 31" x 32"	125 LBS
SMI11	10.9	309	8°C above ambient to 70°C	36.25" x 20" x 26"	42" x 27" x 38"	250 LBS
SMI12	11.7	332	8°C above ambient to 70°C	19.5" x 20" x 26"	25.5" x 27.5" x 76"	250 LBS
SMI31	30.9	874	8°C above ambient to 70°C	32.5" x 26" x 63.5"	39" x 34.5" x 77"	525 LBS
SMI39	38.6	1093	8°C above ambient to 70°C	35.5" x 26" x 76"	41.5" x 34.5" x 90.5"	550 LBS

CO2 INCUBATOR

SCO2W	1.5	42	8°C above ambient to 60°C	15.5" x 15.5" x 10"	21" x 22.5" x 27"	150 LBS
SCO5W	5	141	8°C above ambient to 60°C	20.5" x 19.75" x 21.5"	27.5" x 28.5" x 38"	275 LBS
SCO5A	5	141	8°C above ambient to 60°C	19.25" x 19.75" x 23"	26" x 25.5" x 40.5"	275 LBS
SCO6WE	5.9	166	8°C above ambient to 60°C	20.25" x 19.75" x 25.25"	26" x 26.5" x 40.5"	325 LBS
SCO31	31.4	890	8°C above ambient to 60°C	32.75" x 26" x 63.25"	39.75" x 33.75" x 76.25"	550 LBS
SCO40	39.5	1118	8°C above ambient to 60°C	34.5" x 25.75" x 76"	42.25" x 33.5" x 88"	625 LBS
SCO58	58	1642	8°C above ambient to 60°C	43" x 34.5" x 67.5"	51" x 45" x 81"	1025 LBS

REFRIGERATED INCUBATOR

SRI3	2.4	68	0°C to 45°C	16" x 14" x 22.5"	24" x 21" x 34"	90 LBS
SRI20P	19.3	546	15°C to 30°C	25.5" x 24" x 54.5"	30" x 31.5" x 69.5"	210 LBS