

CERVICAL FLARED EDGE COATED HEADREST

PRODUCT INFORMATION

PRODUCT INFORMATION
CERVICAL FLARED EDGE COATED HEADREST
PPFP-PCAR



The Cervical Flared Edge Coated Headrest supports and positions the head and cervical area during medical imaging procedures. Measuring 16" L x 7.5" W x 4.5" H, its contoured design provides support for the head and neck during examinations.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the Cervical Flared Edge Coated Headrest is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.


SPECIFICATIONS

SHIPPING WEIGHT: 32 oz	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: 17 x 9 x 9 in	POSITIONER TYPE: Neck & Head Positioners
PRODUCT DIMENSIONS: 16" L x 7.5" W x 4.5" H	FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

CERVICAL FLARED EDGE COATED HEADREST

PPFP-PCAR



Date: 9/21/2026		Drawing Title: Cervical Flared Edge Coated Headrest
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAR
		Notes: All dimensions are ±0.1"