

BARIATRIC TORSO FLARED EDGE COATED BOLSTER

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION
BARIATRIC TORSO FLARED EDGE COATED BOLSTER
PPFP-PCAT



The Bariatric Torso Flared Edge Coated Bolster supports and positions the torso during medical imaging procedures. Measuring 18.5" L x 16.5" W x 7" H, its angled design provides a broad positioning surface for bariatric patients during examinations.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the Bariatric Torso Flared Edge Coated Bolster is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.


SPECIFICATIONS

SHIPPING WEIGHT: 64 oz	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: 28 x 16 x 12 in	POSITIONER TYPE: Torso Positioners
PRODUCT DIMENSIONS: 18.5" L x 16.5" W x 7" H	FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

BARIATRIC TORSO FLARED EDGE COATED BOLSTER

PPFP-PCAT



Date: 9/21/2026		Drawing Title: Bariatric Torso Flared Edge Coated Bolster
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAT
		Notes: All dimensions are ±0.1"