


45 DEGREE BARIATRIC SPINAL FLARED EDGE COATED WEDGE

PRODUCT INFORMATION

PRODUCT INFORMATION

45 DEGREE BARIATRIC SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBA



The 45 Degree Bariatric Spinal Flared Edge Coated Wedge supports and positions the body during spinal imaging procedures. Measuring 27.5" L x 10" W x 7" H, its 45-degree wedge design provides angled support for bariatric positioning needs.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the 45 Degree Bariatric Spinal Flared Edge Coated Wedge is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 64 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in

PRODUCT DIMENSIONS: 27.5" L x 10" W x 7" H

MATERIAL: Coated Foam


POSITIONER TYPE: Positioning Wedges

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

45 DEGREE BARIATRIC SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBA



Date: 9/21/2026		Drawing Title: 45 Degree Bariatric Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBA
		Notes: All dimensions are ± 0.1 "