

45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PRODUCT INFORMATION

PRODUCT INFORMATION
45 DEGREE SPINAL FLARED EDGE COATED WEDGE
PPFP-PCBC | PPFP-PCBB | PPFP-PCBD | PPFP-PCEB



The 45 Degree Spinal Flared Edge Coated Wedge supports and positions the body during spinal imaging procedures. Its 45-degree design provides angled support for positioning during examinations. The wedge is available in four sizes:

- **Small:** 13.5" L × 10" W × 7" H
- **Medium:** 20.25" L × 10" W × 7" H
- **Large:** 27.75" L × 10" W × 7" H
- **Extra Large:** 30" L × 10" W × 7" H

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the 45 Degree Spinal Flared Edge Coated Wedge is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.


SPECIFICATIONS

SHIPPING WEIGHT: Multiple options	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: Multiple options	POSITIONER TYPE: Positioning Wedges
PRODUCT SIZES: Small, Medium, Large, XL	FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBC




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBC
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBB




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBB
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBD




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBD
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCEB



Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCEB
		Notes: All dimensions are ± 0.1 "