

FLARED EDGE COATED WEDGE

PRODUCT INFORMATION

PRODUCT INFORMATION
FLARED EDGE COATED WEDGE
PPFP-PCBF | PPFP-PCBG



The Flared Edge Coated Wedge provides angled support for patient positioning during medical imaging procedures. It is available in two angle and size options:

- **21 Degree:** 14.25" L × 12.25" W × 4.5" H
- **27 Degree:** 15.5" L × 12.75" W × 6" H

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.


The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the Flared Edge Coated Wedge is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS


SHIPPING WEIGHT: Multiple options	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: Multiple options	POSITIONER TYPE: Positioning Wedges
AVAILABLE ANGLES: 21 and 27 degrees	FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
FLARED EDGE COATED WEDGE
PPFP-PCBF



Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBF
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION
FLARED EDGE COATED WEDGE
PPFP-PCBG



Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBG
		Notes: All dimensions are ± 0.1 "