

30-60-90 MULTI-ANGLE FLARED EDGE COATED WEDGE

PRODUCT INFORMATION

PRODUCT INFORMATION
30-60-90 MULTI-ANGLE FLARED EDGE COATED WEDGE
PPFP-PCBU



The 30-60-90 Multi-Angle Flared Edge Coated Wedge provides angled support for patient positioning during medical imaging procedures. Measuring 29.5" L x 10" W x 6" H, its multi-angle design provides 30, 60, and 90 degree positioning options within a single wedge.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the 30-60-90 Multi-Angle Flared Edge Coated Wedge is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 64 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in


PRODUCT DIMENSIONS: 29.5" L x 10" W x 6" H

MATERIAL: Coated Foam

POSITIONER TYPE: Positioning Wedges

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
30-60-90 MULTI-ANGLE FLARED EDGE COATED WEDGE
PPFP-PCBU



Date: 9/21/2026		Drawing Title: 30-60-90 Multi-Angle Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBU
		Notes: All dimensions are ± 0.1 "