


HEAD & NECK FLARED EDGE COATED SUPPORT

PRODUCT INFORMATION

PRODUCT INFORMATION
HEAD & NECK FLARED EDGE COATED SUPPORT
PPFP-PCAL



The Head & Neck Flared Edge Coated Support provides support and stabilization for the head and neck during medical imaging procedures. Measuring 16" L x 14.5" W x 6.5" H, its contoured design helps maintain positioning while supporting the head and neck.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design uses flared edges to significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size, they do not affect imaging or positioning.

Finally, the Head & Neck Flared Edge Coated Support is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS


SHIPPING WEIGHT: 48 oz
SHIPPING DIMENSIONS: 26 x 11 x 14 in
PRODUCT DIMENSIONS: 16" L x 14.5" W x 6.5" H

MATERIAL: Coated Foam
POSITIONER TYPE: Neck & Head Positioners
FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

HEAD & NECK FLARED EDGE COATED SUPPORT

PPFP-PCAL



Date: 9/21/2026		Drawing Title: Head & Neck Flared Edge Coated Support
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAL
		Notes: All dimensions are ± 0.1 "