

COATED RECTANGLE POSITIONER

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION

COATED RECTANGLE POSITIONER

PPFP-PCBR | PPFP-PCBS | PPFP-PCBT | PPFP-PCBQ | PPFP-PCEM



The Coated Rectangle Positioner provides a flat, stable surface for patient positioning during medical imaging procedures. It is available in five size options to accommodate different positioning needs:

- **Extra Small:** 6" L × 4" W × 0.5" H
- **Small:** 6" L × 4" W × 1" H
- **Medium:** 6" L × 4" W × 1.5" H
- **Large:** 6" L × 4" W × 2" H
- **Extra Large:** 24" L × 15" W × 2.5" H

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.


These rectangle positioners use a standard non-flared design. As a result, the coated edges may produce visible artifacts during imaging.

Finally, these Coated Rectangle Positioners are MRI Safe and latex-free. This makes them suitable for use across a variety of medical imaging environments.


SPECIFICATIONS

SHIPPING WEIGHT: Multiple options	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: Multiple options	POSITIONER TYPE: Positioning Rectangles
PRODUCT SIZES: XS, Small, Medium, Large, XL	FEATURES: MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
COATED RECTANGLE POSITIONER
PPFP-PCBQ




Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBQ
		Notes: All dimensions are $\pm 0.1"$




Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBR
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION
COATED RECTANGLE POSITIONER
PPFP-PCBS




Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBS
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION
COATED RECTANGLE POSITIONER
PPFP-PCBT



Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBT
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION
COATED RECTANGLE POSITIONER
PPFP-PCEM



Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCEM
		Notes: All dimensions are ± 0.1 "