


ENDO-ULTRASOUND COATED WEDGE

PRODUCT INFORMATION

PRODUCT INFORMATION
ENDO-ULTRASOUND COATED WEDGE
PPFP-PCUZ



The Endo-Ultrasound Coated Wedge provides support for patient positioning during endoscopic ultrasound procedures. Measuring 20" L x 20" W x 6" H, its contoured design includes a recessed area to help accommodate and support the patient during positioning.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

This positioner uses a standard non-flared design. As a result, the coated edges may produce visible artifacts during imaging.

Finally, the Endo-Ultrasound Coated Wedge is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS


SHIPPING WEIGHT: 96 oz
SHIPPING DIMENSIONS: 26 x 24 x 9 in
PRODUCT DIMENSIONS: 20" L x 20" W x 6" H

MATERIAL: Coated Foam
POSITIONER TYPE: Positioning Wedges
FEATURES: MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

ENDO-ULTRASOUND COATED WEDGE

PPFP-PCUZ



Date: 9/21/2026		Drawing Title: Endo-Ultrasound Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCUZ
		Notes: All dimensions are $\pm 0.1"$