


CT COATED HEADREST

PRODUCT INFORMATION

PRODUCT INFORMATION
CT COATED HEADREST
PPFP-PCBW



The CT Coated Headrest provides support and stabilization for the head during CT imaging procedures. Measuring 12.5" L × 11" W × 10" H, its angled design features a contoured recess to help support the head during positioning.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the headrest is radiolucent, with low attenuation and minimal shadowing during imaging.

This headrest uses a standard non-flared design. As a result, it is not guaranteed to be artifact-free, and the coated edges may produce visible artifacts during imaging.

Finally, the CT Coated Headrest is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.


SPECIFICATIONS

SHIPPING WEIGHT: 48 oz	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: 26 × 14 × 11 in	POSITIONER TYPE: Neck & Head Positioners
PRODUCT DIMENSIONS: 12.5" L × 11" W × 10" H	FEATURES: MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

CT COATED HEADREST

PPFP-PCBW



Date: 9/21/2026		Drawing Title: CT Coated Headrest
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBW
		Notes: All dimensions are ±0.1"