


CHILD MYELOGRAM COATED BLOCK

PRODUCT INFORMATION

PRODUCT INFORMATION
CHILD MYELOGRAM COATED BLOCK
PPFP-PCAX



The Child Myelogram Coated Block provides support and stabilization for pediatric patient positioning during medical imaging procedures. Measuring 14" L x 8" W x 10" H, its contoured design helps support the patient and maintain positioning during imaging.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.


This positioner uses a standard non-flared design. As a result, the coated edges may produce visible artifacts during imaging.

Finally, the Child Myelogram Coated Block is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 48 oz	MATERIAL: Coated Foam
SHIPPING DIMENSIONS: 26 x 11 x 14 in	POSITIONER TYPE: Torso Positioners, Pediatric
PRODUCT DIMENSIONS: 14" L x 8" W x 10" H	FEATURES: MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
CHILD MYELOGRAM COATED BLOCK
PPFP-PCAX



Date: 9/21/2026		Drawing Title: Child Myelogram Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAX
		Notes: All dimensions are ± 0.1 "