

DECUBITUS ABDOMINAL FLARED EDGE COATED PAD

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION
DECUBITUS ABDOMINAL FLARED EDGE COATED PAD
PPFP-PCBP



The Decubitus Abdominal Flared Edge Coated Pad provides support and positioning for the abdominal area during medical imaging procedures. It measures 24" L x 18" W at the top and 27.5" L x 21.5" W x 4" H at the base. Its large, flat surface provides support for decubitus positioning during imaging.

For added durability, a medical-grade, fluid-resistant coating covers the foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioner is radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. While the flared sides increase the overall size of the positioner, they do not affect imaging or positioning.

Finally, the Decubitus Abdominal Flared Edge Coated Pad is MRI Safe and latex-free. This makes it suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 112 oz

SHIPPING DIMENSIONS: 47 x 24 x 9 in


PRODUCT DIMENSIONS:
Top: 24" L x 18" W / Base: 27.5" L x 21.5" W x 4" H

MATERIAL: Coated Foam

POSITIONER TYPE: Torso Positioners

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
DECUBITUS ABDOMINAL FLARED EDGE COATED PAD
PPFP-PCBP



Date: 9/21/2026		Drawing Title: Decubitus Abdominal Flared Edge Coated Pad
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBP
		Notes: All dimensions are $\pm 0.1"$