


CLINIC BUNDLE A FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION
CLINIC BUNDLE A FLARED EDGE COATED POSITIONERS
PPFP-KIT2



The Clinic Bundle A Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for clinical imaging environments, it includes blocks and wedges for a range of positioning needs.

The Clinic Bundle A includes **one of each** of the following:

- PPFP-PCAA – Adult Finger Oblique Flared Edge Coated Block
- PPFP-PCAK – Circular Head Flared Edge Coated Block
- PPFP-PCBB – Medium 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBI – 15 Degree Flared Edge Coated Wedge Set

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in Clinic Bundle A are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 64 oz


SHIPPING DIMENSIONS: 26 x 14 x 11 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION
ADULT FINGER OBLIQUE FLARED EDGE COATED BLOCK
PPFP-PCAA




Date: 9/21/2026		Drawing Title: Adult Finger Oblique Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION


CIRCULAR HEAD FLARED EDGE COATED BLOCK

PPFP-PCAK



Date: 9/21/2026		Drawing Title: Circular Head Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAK
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION
MEDIUM 45 DEGREE SPINAL FLARED EDGE COATED WEDGE
PPFP-PCBB





Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBB
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

15 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBI



Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBI
		Notes: All dimensions are ± 0.1 "