

GENERAL BUNDLE C FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION

GENERAL BUNDLE C FLARED EDGE COATED POSITIONERS

PPFP-KIT7



The General Bundle C Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for general imaging environments, it includes blocks, head supports, a bolster, and a rectangle positioner for a range of positioning needs.

The General Bundle C includes **one of each** of the following:

- PPFP-PCAA – Adult Finger Oblique Flared Edge Coated Block
- PPFP-PCAH – Knee Arthrogram Flared Edge Coated Block
- PPFP-PCAM – Adult Flared Edge Coated Head Holder
- PPFP-PCAR – Cervical Flared Edge Coated Head Rest
- PPFP-PCAY – Angular Flared Edge Coated Bolster
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in General Bundle C are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 160 oz

SHIPPING DIMENSIONS: 15 × 15 × 41 in


MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION
ADULT FINGER OBLIQUE FLARED EDGE COATED BLOCK
PPFP-PCAA




Date: 9/21/2026		Drawing Title: Adult Finger Oblique Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

KNEE ARTHROGRAM FLARED EDGE COATED BLOCK

PPFP-PCAH




Date: 9/21/2026		Drawing Title: Knee Arthrogram Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAH
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

ADULT FLARED EDGE COATED HEAD HOLDER

PPFP-PCAM




Date: 9/21/2026		Drawing Title: Adult Flared Edge Coated Head Holder
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAM
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

CERVICAL FLARED EDGE COATED HEAD REST

PPFP-PCAR




Date: 9/21/2026		Drawing Title: Cervical Flared Edge Coated Headrest
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAR
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

ANGULAR FLARED EDGE COATED BOLSTER

PPFP-PCAY




Date: 9/21/2026		Drawing Title: Angular Flared Edge Coated Bolster
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAY
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL



Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are $\pm 0.1"$