

GENERAL BUNDLE E FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION
GENERAL BUNDLE E FLARED EDGE COATED POSITIONERS
PPFP-KIT9



The General Bundle E Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for general imaging environments, it includes blocks, a head holder, wedges, a wedge set, and rectangle positioners for a range of positioning needs.

The General Bundle E includes **one of each** of the following:

- PPFP-PCAA – Adult Finger Oblique Flared Edge Coated Block
- PPFP-PCAK – Circular Head Flared Edge Coated Block
- PPFP-PCAM – Adult Flared Edge Coated Head Holder
- PPFP-PCBB – Medium 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBG – 21 Degree Flared Edge Coated Wedge
- PPFP-PCBI – 15 Degree Flared Edge Coated Wedge Set
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner
- PPFP-PCBN – Large 3" Flared Edge Coated Rectangle Positioner
- PPFP-PCBU – 30-60-90 Multi-Angle Flared Edge Coated Wedge

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in General Bundle E are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS


SHIPPING WEIGHT: 176 oz
SHIPPING DIMENSIONS: 14 x 14 x 40 in
MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets
FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

ADULT FINGER OBLIQUE FLARED EDGE COATED BLOCK

PPFP-PCAA




Date: 9/21/2026		Drawing Title: Adult Finger Oblique Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

CIRCULAR HEAD FLARED EDGE COATED BLOCK

PPFP-PCAK




Date: 9/21/2026		Drawing Title: Circular Head Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAK
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION


ADULT FLARED EDGE COATED HEAD HOLDER

PPFP-PCAM



Date: 9/21/2026		Drawing Title: Adult Flared Edge Coated Head Holder
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAM
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION
MEDIUM 45 DEGREE SPINAL FLARED EDGE COATED WEDGE
PPFP-PCBB




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBB
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

21 DEGREE FLARED EDGE COATED WEDGE

PPFP-PCBG





Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBG
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

15 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBI




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBI
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION


SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL




Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION
LARGE 3" FLARED EDGE COATED RECTANGLE POSITIONER
PPFP-PCBN



Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBN
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION
30-60-90 MULTI-ANGLE FLARED EDGE COATED WEDGE
PPFP-PCBU



Date: 9/21/2026		Drawing Title: 30-60-90 Multi-Angle Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBU
		Notes: All dimensions are ± 0.1 "