

GENERAL BUNDLE G FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION

GENERAL BUNDLE G FLARED EDGE COATED POSITIONERS

PPFP-KIT10



The General Bundle G Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for general imaging environments, it includes blocks, wedges, wedge sets, and a rectangle positioner for a range of positioning needs.

The General Bundle G includes **one of each** of the following:

- PPFP-PCAA – Adult Finger Oblique Flared Edge Coated Block
- PPFP-PCAK – Circular Head Flared Edge Coated Block
- PPFP-PCBD – Large 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBE – 34 Degree Flared Edge Coated Wedge Set
- PPFP-PCBI – 15 Degree Flared Edge Coated Wedge Set
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in General Bundle G are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 80 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

ADULT FINGER OBLIQUE FLARED EDGE COATED BLOCK

PPFP-PCAA




Date: 9/21/2026		Drawing Title: Adult Finger Oblique Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

CIRCULAR HEAD FLARED EDGE COATED BLOCK

PPFP-PCAK




Date: 9/21/2026		Drawing Title: Circular Head Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAK
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

LARGE 45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBD





Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBD
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

34 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBE





Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBE
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

15 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBI




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBI
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL



Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are ± 0.1 "