

GENERAL BUNDLE 1 FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION
GENERAL BUNDLE 1 FLARED EDGE COATED POSITIONERS
PPFP-KIT12



The General Bundle 1 Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for general imaging environments, it includes head blocks, wedges, and rectangle positioners for a range of positioning needs.

The General Bundle 1 includes **one of each** of the following, except where otherwise specified:

- PPFP-PCAK – Circular Head Flared Edge Coated Block (Qty. 2)
- PPFP-PCBD – Large 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBI – 15 Degree Flared Edge Coated Wedge Set
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner
- PPFP-PCBT – Large Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. The flared edge positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design featured on the flared edge positioners helps significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning. The Coated Rectangle Positioner has a traditional non-flared design and is not radiolucent.

Finally, the positioners included in General Bundle 1 are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 80 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

CIRCULAR HEAD FLARED EDGE COATED BLOCK

PPFP-PCAK




Date: 9/21/2026		Drawing Title: Circular Head Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAK
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

LARGE 45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBD





Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBD
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

15 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBI




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBI
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL




Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

LARGE COATED RECTANGLE POSITIONER

PPFP-PCBT



Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBT
		Notes: All dimensions are ±0.1"