


ORTHOPEDIC BUNDLE A FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION

ORTHOPEDIC BUNDLE A FLARED EDGE COATED POSITIONERS

PPFP-KIT13



The Orthopedic Bundle A Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for orthopedic imaging environments, it includes a block, head supports, a wedge, and a rectangle positioner for a range of positioning needs.

The Orthopedic Bundle A includes **one of each** of the following:

- PPFP-PCAA – Adult Finger Oblique Flared Edge Coated Block
- PPFP-PCAM – Adult Flared Edge Coated Head Holder
- PPFP-PCAR – Cervical Flared Edge Coated Head Rest
- PPFP-PCBD – Large 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in Orthopedic Bundle A are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 96 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

ADULT FINGER OBLIQUE FLARED EDGE COATED BLOCK

PPFP-PCAA




Date: 9/21/2026		Drawing Title: Adult Finger Oblique Flared Edge Coated Block
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

ADULT FLARED EDGE COATED HEAD HOLDER

PPFP-PCAM




Date: 9/21/2026		Drawing Title: Adult Flared Edge Coated Head Holder
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAM
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

CERVICAL FLARED EDGE COATED HEAD REST

PPFP-PCAR




Date: 9/21/2026		Drawing Title: Cervical Flared Edge Coated Headrest
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAR
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

LARGE 45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBD




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBD
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL



Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are ± 0.1 "