

ORTHOPEDIC BUNDLE B FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION

ORTHOPEDIC BUNDLE B FLARED EDGE COATED POSITIONERS

PPFP-KIT14



The Orthopedic Bundle B Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for orthopedic imaging environments, it includes head and neck supports, wedges, a wedge set, and a rectangle positioner for a range of positioning needs.

The Orthopedic Bundle B includes **one of each** of the following:

- PPFP-PCAL – Head & Neck Flared Edge Coated Support
- PPFP-PCAN – Child Flared Edge Coated Head Holder
- PPFP-PCAR – Cervical Flared Edge Coated Head Rest
- PPFP-PCBC – Small 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBG – 21 Degree Flared Edge Coated Wedge
- PPFP-PCBJ – 12 Degree Flared Edge Coated Wedge Set
- PPFP-PCBL – Small 3" Flared Edge Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in Orthopedic Bundle B are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 112 oz

SHIPPING DIMENSIONS: 26 x 24 x 9 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

HEAD & NECK FLARED EDGE COATED SUPPORT

PPFP-PCAL




Date: 9/21/2026		Drawing Title: Head & Neck Flared Edge Coated Support
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAL
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

CHILD FLARED EDGE COATED HEAD HOLDER

PPFP-PCAN




Date: 9/21/2026		Drawing Title: Child Flared Edge Coated Head Holder
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAN
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

CERVICAL FLARED EDGE COATED HEAD REST

PPFP-PCAR




Date: 9/21/2026		Drawing Title: Cervical Flared Edge Coated Headrest
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAR
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

SMALL 45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBC




Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBC
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

21 DEGREE FLARED EDGE COATED WEDGE

PPFP-PCBG





Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBG
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

12 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBJ




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBJ
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

SMALL 3" FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBL



Date: 9/21/2026		Drawing Title: 3" Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBL
		Notes: All dimensions are ±0.1"