

WEDGE BUNDLE B FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION

WEDGE BUNDLE B FLARED EDGE COATED POSITIONERS

PPFP-KIT18



The Wedge Bundle B Flared Edge Coated Positioners is a kit that brings together a selection of foam wedges for patient support and positioning during medical imaging procedures. Designed to provide a variety of positioning angles, it includes wedges and a wedge set for a range of positioning needs.

The Wedge Bundle B includes **one of each** of the following:

- PPFP-PCBD – Large 45 Degree Spinal Flared Edge Coated Wedge
- PPFP-PCBE – 34 Degree Flared Edge Coated Wedge Set
- PPFP-PCBF – 27 Degree Flared Edge Coated Wedge
- PPFP-PCBG – 21 Degree Flared Edge Coated Wedge

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in Wedge Bundle B are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 80 oz

SHIPPING DIMENSIONS: 28 x 16 x 12 in

MATERIAL: Coated Foam


POSITIONER TYPE: X-Ray Positioning Sets

FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

LARGE 45 DEGREE SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBD





Date: 9/21/2026		Drawing Title: 45 Degree Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBD
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION

34 DEGREE FLARED EDGE COATED WEDGE SET

PPFP-PCBE




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge Set
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBE
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

27 DEGREE FLARED EDGE COATED WEDGE

PPFP-PCBF




Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBF
		Notes: All dimensions are ±0.1"

PRODUCT INFORMATION

21 DEGREE FLARED EDGE COATED WEDGE

PPFP-PCBG



Date: 9/21/2026		Drawing Title: Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBG
		Notes: All dimensions are ±0.1"