

BARIATRIC BUNDLE FLARED EDGE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION
BARIATRIC BUNDLE FLARED EDGE COATED POSITIONERS
PPFP-KIT19



The Bariatric Bundle Flared Edge Coated Positioners is a kit that brings together a selection of foam positioners for patient support and positioning during medical imaging procedures. Designed for bariatric positioning needs, it includes a bolster, wedge, rectangle positioner, and abdominal pad for a range of positioning needs.

The Bariatric Bundle includes **one of each** of the following:

- PPFP-PCAT – Bariatric Torso Flared Edge Coated Bolster
- PPFP-PCBA – 45 Degree Bariatric Spinal Flared Edge Coated Wedge
- PPFP-PCBM – Bariatric Flared Edge Coated Rectangle Positioner
- PPFP-PCBP – Decubitus Abdominal Flared Edge Coated Pad

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The patented N-VISI-RAY design features flared edges that help significantly reduce artifacts during X-ray imaging. Although the flared sides increase the overall size of the positioners, they do not affect imaging or positioning.

Finally, the positioners included in the Bariatric Bundle are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 176 oz

SHIPPING DIMENSIONS: 47 x 24 x 9 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: Flared Edge, MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

BARIATRIC TORSO FLARED EDGE COATED BOLSTER

PPFP-PCAT




Date: 9/21/2026		Drawing Title: Bariatric Torso Flared Edge Coated Bolster
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCAT
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION

45 DEGREE BARIATRIC SPINAL FLARED EDGE COATED WEDGE

PPFP-PCBA




Date: 9/21/2026		Drawing Title: 45 Degree Bariatric Spinal Flared Edge Coated Wedge
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBA
		Notes: All dimensions are ± 0.1 "

PRODUCT INFORMATION


BARIATRIC FLARED EDGE COATED RECTANGLE POSITIONER

PPFP-PCBM



Date: 9/21/2026		Drawing Title: Bariatric Flared Edge Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBM
		Notes: All dimensions are $\pm 0.1"$

PRODUCT INFORMATION
DECUBITUS ABDOMINAL FLARED EDGE COATED PAD
PPFP-PCBP



Date: 9/21/2026		Drawing Title: Decubitus Abdominal Flared Edge Coated Pad
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBP
		Notes: All dimensions are $\pm 0.1"$