

SCANNING SUPPORT BUNDLE COATED POSITIONERS

PRODUCT INFORMATION

PRODUCT INFORMATION

SCANNING SUPPORT BUNDLE COATED POSITIONERS

PPFP-KIT21



The Scanning Support Bundle Coated Positioners is a kit that brings together foam rectangle positioners for support during ultrasound scanning. Designed for general and vascular sonographers, the positioners can be used together to support the scanning arm, helping reduce fatigue and improve postural alignment during procedures.

The Scanning Support Bundle includes **ten** of the following:

- PPFP-PCBT – Large Coated Rectangle Positioner

For added durability, a medical-grade, fluid-resistant coating covers each foam positioner. The coating is easy to clean with medical-grade disinfectant wipes or sprays. It also helps protect the foam interior during repeated clinical use. In addition, the positioners are radiolucent, with low attenuation and minimal shadowing during imaging.

The positioners feature traditional non-flared edges. As a result, the edges may produce artifacts during imaging, particularly where the coating is applied.

Finally, the positioners included in the Scanning Support Bundle are MRI Safe and latex-free. This makes the kit suitable for use across a variety of medical imaging environments.

SPECIFICATIONS

SHIPPING WEIGHT: 32 oz

SHIPPING DIMENSIONS: 17 × 9 × 9 in

MATERIAL: Coated Foam

POSITIONER TYPE: X-Ray Positioning Sets


FEATURES: MRI Safe, Latex-Free, Radiolucent, Low Attenuation, Easy to Clean

PRODUCT INFORMATION

PRODUCT INFORMATION

LARGE COATED RECTANGLE POSITIONER

PPFP-PCBT



Date: 9/21/2026		Drawing Title: Coated Rectangle Positioner
Scale: 1:1	Unit: Inches	Drawing No: PPFP-PCBT
		Notes: All dimensions are ±0.1"