

CANDo® MGE® PLUS exerciser



- CanDo® MGE® PLUS is made of woven polyester with latex-free elastic
- Easy-to-use grips allow for unlimited exercise options
- Each exerciser contains both large and small grips. Use small grips as anchors or for finger exercises
- Wrap around hands and feet to perform high-impact upper and lower body exercises
- Woven exerciser contains no latex, no scent, and no powder coating
- Can be washed and dried
- Easily change resistance by moving 1 grip (loop) up or down the exerciser, or using a different resistance
- Progressively exercise with 8 different color-coded resistance levels available in the Preferred Colors® sequence
- Available in 9-grip (6 ft. exercisers) and 135-grip (30 yd. dispenser rolls) equivalent to (15) 6 ft. exercisers




		9 grips, 6 foot exerciser	135 grips, 30 yard roll, 15 exercisers
<div></div> tan	xx-light	10-6590	10-6760
<div></div> yellow	x-light	10-6591	10-6761
<div></div> red	light	10-6592	10-6762
<div></div> green	medium	10-6593	10-6763
<div></div> blue	heavy	10-6594	10-6764
<div></div> black	x-heavy	10-6595	10-6765
<div></div> silver	xx-heavy	10-6596	10-6766
<div></div> gold	xxx-heavy	10-6597	10-6767
<div></div> yel red grn blu blk	set of 5 (1 ea)	10-6598	10-6768



use for upper body exercises



use for core exercises



use for lower body exercises

