

Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Mini Refill, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10-6
Pack Qty	384	Centrifugation	NA
Case Qty	3840	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Mini Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-250-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Mini Refill

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	384	Centrifugation	NA
Case Qty	3840	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Mini Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-800-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in Resealable Bags

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1000	Centrifugation	NA
Case Qty	10000	Centrifugation Speed (RCF)	
Package Type	Resealable Bags	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-260-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Refill

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	768	Centrifugation	NA
Case Qty	7680	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-265-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Refills, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	768	Centrifugation	NA
Case Qty	7680	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-960-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in 96 Racks

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	96 Rack	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1015-965-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in 96 Racks, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	96 Rack	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 20-Oct-2026 . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023