


# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-250-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Mini Refill

|               |                      |                            |                 |
|---------------|----------------------|----------------------------|-----------------|
| Fluid Measure | μL                   | Freezable                  | NA              |
| Fluid Volume  | 20                   | Freeze Temp                |                 |
| Height (mm)   | 48.3                 | Autoclavable               | 121°C 15min     |
| Width (mm)    |                      | Sterile                    | No              |
| Length (mm)   |                      | Sterile Processing         | NA              |
| Pack Qty      | 480                  | Centrifugation             | NA              |
| Case Qty      | 4800                 | Centrifugation Speed (RCF) |                 |
| Package Type  | Eclipse™ Mini Refill | Shelf Life (Months)        | 72              |
| Material      | Polypropylene        | Certification              | ISO 9001        |
| Cap Material  |                      | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


### Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-255-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Mini Refill, Sterile

|               |                      |                            |                                                                                 |
|---------------|----------------------|----------------------------|---------------------------------------------------------------------------------|
| Fluid Measure | μL                   | Freezable                  | NA                                                                              |
| Fluid Volume  | 20                   | Freeze Temp                |                                                                                 |
| Height (mm)   | 48.3                 | Autoclavable               | 121°C 15min                                                                     |
| Width (mm)    |                      | Sterile                    | Yes                                                                             |
| Length (mm)   |                      | Sterile Processing         | Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup> |
| Pack Qty      | 480                  | Centrifugation             | NA                                                                              |
| Case Qty      | 4800                 | Centrifugation Speed (RCF) |                                                                                 |
| Package Type  | Eclipse™ Mini Refill | Shelf Life (Months)        | 60                                                                              |
| Material      | Polypropylene        | Certification              | ISO 9001                                                                        |
| Cap Material  |                      | Regulatory                 | U.S. FDA 21 CFR                                                                 |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-800-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Resealable Bags

|               |                 |                            |                 |
|---------------|-----------------|----------------------------|-----------------|
| Fluid Measure | μL              | Freezable                  | NA              |
| Fluid Volume  | 20              | Freeze Temp                |                 |
| Height (mm)   | 48.3            | Autoclavable               | 121°C 15min     |
| Width (mm)    |                 | Sterile                    | No              |
| Length (mm)   |                 | Sterile Processing         | NA              |
| Pack Qty      | 1000            | Centrifugation             | NA              |
| Case Qty      | 10000           | Centrifugation Speed (RCF) |                 |
| Package Type  | Resealable Bags | Shelf Life (Months)        | 72              |
| Material      | Polypropylene   | Certification              | ISO 9001        |
| Cap Material  |                 | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-260-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Refill

|               |                 |                            |                 |
|---------------|-----------------|----------------------------|-----------------|
| Fluid Measure | μL              | Freezable                  | NA              |
| Fluid Volume  | 20              | Freeze Temp                |                 |
| Height (mm)   | 48.3            | Autoclavable               | 121°C 15min     |
| Width (mm)    |                 | Sterile                    | No              |
| Length (mm)   |                 | Sterile Processing         | NA              |
| Pack Qty      | 1056            | Centrifugation             | NA              |
| Case Qty      | 10560           | Centrifugation Speed (RCF) |                 |
| Package Type  | Eclipse™ Refill | Shelf Life (Months)        | 72              |
| Material      | Polypropylene   | Certification              | ISO 9001        |
| Cap Material  |                 | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


### Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-290-020-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Pagoda® Refill

|               |                |                            |                 |
|---------------|----------------|----------------------------|-----------------|
| Fluid Measure | μL             | Freezable                  | NA              |
| Fluid Volume  | 20             | Freeze Temp                |                 |
| Height (mm)   | 48.3           | Autoclavable               | 121°C 15min     |
| Width (mm)    |                | Sterile                    | No              |
| Length (mm)   |                | Sterile Processing         | NA              |
| Pack Qty      | 1056           | Centrifugation             | NA              |
| Case Qty      | 10560          | Centrifugation Speed (RCF) |                 |
| Package Type  | Pagoda® Refill | Shelf Life (Months)        | 72              |
| Material      | Polypropylene  | Certification              | ISO 9001        |
| Cap Material  |                | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-700-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Stack Racks

|               |                |                            |                 |
|---------------|----------------|----------------------------|-----------------|
| Fluid Measure | μL             | Freezable                  | NA              |
| Fluid Volume  | 20             | Freeze Temp                |                 |
| Height (mm)   | 48.3           | Autoclavable               | 121°C 15min     |
| Width (mm)    |                | Sterile                    | No              |
| Length (mm)   |                | Sterile Processing         | NA              |
| Pack Qty      | 960            | Centrifugation             | NA              |
| Case Qty      | 9600           | Centrifugation Speed (RCF) |                 |
| Package Type  | 192 Stack Rack | Shelf Life (Months)        | 72              |
| Material      | Polypropylene  | Certification              | ISO 9001        |
| Cap Material  |                | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-265-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in Eclipse™ Refills, Sterile

|               |                 |                            |                                                                                 |
|---------------|-----------------|----------------------------|---------------------------------------------------------------------------------|
| Fluid Measure | μL              | Freezable                  | NA                                                                              |
| Fluid Volume  | 20              | Freeze Temp                |                                                                                 |
| Height (mm)   | 48.3            | Autoclavable               | 121°C 15min                                                                     |
| Width (mm)    |                 | Sterile                    | Yes                                                                             |
| Length (mm)   |                 | Sterile Processing         | Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup> |
| Pack Qty      | 1056            | Centrifugation             | NA                                                                              |
| Case Qty      | 10560           | Centrifugation Speed (RCF) |                                                                                 |
| Package Type  | Eclipse™ Refill | Shelf Life (Months)        | 60                                                                              |
| Material      | Polypropylene   | Certification              | ISO 9001                                                                        |
| Cap Material  |                 | Regulatory                 | U.S. FDA 21 CFR                                                                 |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


### Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-960-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in 96 Racks

|               |               |                            |                 |
|---------------|---------------|----------------------------|-----------------|
| Fluid Measure | μL            | Freezable                  | NA              |
| Fluid Volume  | 20            | Freeze Temp                |                 |
| Height (mm)   | 48.3          | Autoclavable               | 121°C 15min     |
| Width (mm)    |               | Sterile                    | No              |
| Length (mm)   |               | Sterile Processing         | NA              |
| Pack Qty      | 960           | Centrifugation             | NA              |
| Case Qty      | 9600          | Centrifugation Speed (RCF) |                 |
| Package Type  | 96 Rack       | Shelf Life (Months)        | 72              |
| Material      | Polypropylene | Certification              | ISO 9001        |
| Cap Material  |               | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-790-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Racks

|               |               |                            |                 |
|---------------|---------------|----------------------------|-----------------|
| Fluid Measure | μL            | Freezable                  | NA              |
| Fluid Volume  | 20            | Freeze Temp                |                 |
| Height (mm)   | 48.3          | Autoclavable               | 121°C 15min     |
| Width (mm)    |               | Sterile                    | No              |
| Length (mm)   |               | Sterile Processing         | NA              |
| Pack Qty      | 960           | Centrifugation             | NA              |
| Case Qty      | 9600          | Centrifugation Speed (RCF) |                 |
| Package Type  | 192 Rack      | Shelf Life (Months)        | 72              |
| Material      | Polypropylene | Certification              | ISO 9001        |
| Cap Material  |               | Regulatory                 | U.S. FDA 21 CFR |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-965-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in 96 Racks, Sterile

|               |               |                            |                                                                                 |
|---------------|---------------|----------------------------|---------------------------------------------------------------------------------|
| Fluid Measure | μL            | Freezable                  | NA                                                                              |
| Fluid Volume  | 20            | Freeze Temp                |                                                                                 |
| Height (mm)   | 48.3          | Autoclavable               | 121°C 15min                                                                     |
| Width (mm)    |               | Sterile                    | Yes                                                                             |
| Length (mm)   |               | Sterile Processing         | Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup> |
| Pack Qty      | 960           | Centrifugation             | NA                                                                              |
| Case Qty      | 9600          | Centrifugation Speed (RCF) |                                                                                 |
| Package Type  | 96 Rack       | Shelf Life (Months)        | 60                                                                              |
| Material      | Polypropylene | Certification              | ISO 9001                                                                        |
| Cap Material  |               | Regulatory                 | U.S. FDA 21 CFR                                                                 |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


### Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

# Technical Product Data

Technical Data Sheet

Labcon Product Number 1011-795-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## Eclipse™ 20 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Racks, Sterile

|               |               |                            |                                                                                 |
|---------------|---------------|----------------------------|---------------------------------------------------------------------------------|
| Fluid Measure | μL            | Freezable                  | NA                                                                              |
| Fluid Volume  | 20            | Freeze Temp                |                                                                                 |
| Height (mm)   | 48.3          | Autoclavable               | 121°C 15min                                                                     |
| Width (mm)    |               | Sterile                    | Yes                                                                             |
| Length (mm)   |               | Sterile Processing         | Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup> |
| Pack Qty      | 960           | Centrifugation             | NA                                                                              |
| Case Qty      | 9600          | Centrifugation Speed (RCF) |                                                                                 |
| Package Type  | 192 Rack      | Shelf Life (Months)        | 60                                                                              |
| Material      | Polypropylene | Certification              | ISO 9001                                                                        |
| Cap Material  |               | Regulatory                 | U.S. FDA 21 CFR                                                                 |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

**Additional Testing if applicable:**

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 20-Oct-2026 . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023