

Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.


Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Stack Racks

Fluid Measure	μL	Freezable	NA
Fluid Volume	300	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	192 Stack Rack	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.


Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Racks

Fluid Measure	μL	Freezable	NA
Fluid Volume	200	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	192 Rack	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 300 uL Pipet Tips for Rainin® LTS Pipettors, in 192 Racks, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	200	Freeze Temp	
Height (mm)	55.6	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10-6
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	192 Rack	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 23-Oct-2026 . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023