


Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-900-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.


Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Bulk

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	38.1	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	30000	Centrifugation	NA
Case Qty	30000	Centrifugation Speed (RCF)	
Package Type	Bulk	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Eclipse™ Mini Refills, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10-6
Pack Qty	576	Centrifugation	NA
Case Qty	5760	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Mini Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-250-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Eclipse™ Mini Refills

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	576	Centrifugation	NA
Case Qty	5760	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Mini Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-800-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Resealable Bags

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	38.1	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1000	Centrifugation	NA
Case Qty	10000	Centrifugation Speed (RCF)	
Package Type	Resealable Bags	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-290-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Pagoda®

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	Pagoda® Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-260-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Eclipse™ Refills

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1248	Centrifugation	NA
Case Qty	12480	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-960-038-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in 96 Racks

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	96 Rack	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-290-020-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Pagoda® Express Refills

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	Pagoda® Express Refill	Shelf Life (Months)	72
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-295-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Pagoda® Refills, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	Pagoda® Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-965-038-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in 96 Racks, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	96 Rack	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-265-000-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Eclipse™ Refills, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	1248	Centrifugation	NA
Case Qty	12480	Centrifugation Speed (RCF)	
Package Type	Eclipse™ Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.


Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

Technical Product Data

Technical Data Sheet

Labcon Product Number 1038-295-020-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ 10 uL Graduated Pipet Tips with UltraFine™ Point, in Pagoda® Express Refills, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	10	Freeze Temp	
Height (mm)	31.4	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	1152	Centrifugation	NA
Case Qty	11520	Centrifugation Speed (RCF)	
Package Type	Pagoda® Express Refill	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 20-Oct-2026 . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023