

Item	Pt #	Pt #	Description	Size	QTY.
1	7822SS	137822000	Nested Table 7822SS	48"L x 24"W x 42"H	1
2	7823SS	137823000	Nested Table 7823SS	44"L x 22"W x 40"H	1
3	7824SS	137824000	Nested Table 7824SS	40"L x 22"W x 38"H	1
4	7825SS	137825000	Nested Table 7825SS	36"L x 18"W x 36"H	1
5	7826SS	137826000	Nested Table 7826SS	32"L x 16"W x 34"H	1
6	7827SS	137827000	Nested Table 7827SS	28"L x 14"W x 32"H	1

NESTED INSTRUMENT TABLE
 are a space saving alternative with a unique design that requires only 8sq.of floor space but yields 31 sq.ft. of usable workspace. They can be order individually or as complete set.

PO#: _____ DATE: _____

APPROVAL: _____

TITLE: _____



PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF BLICKMAN HEALTH INDUSTRIES. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF BLICKMAN HEALTH INDUSTRIES IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:
 DIMENSIONS ARE IN INCHES
 BEND RADIUS = MTL. THK.
 TOLERANCES:
 FRACTION: $\pm 1/64"$
 DECIMAL: $\pm .015"$

	NAME	DATE
DRAWN	AP	04/12/2006
ENG APPR.		
MATERIAL	N/A	
FINISH	N/A	
APPRX WEIGHT = 177.423 lbs		
DO NOT SCALE DRAWING		

BLICKMAN INC.
 CLIFTON, N.J.

TITLE:

Nested Table Set
 7821SS

SIZE
A

DWG. NO.
 0137821000





REV
 PC


SCALE: 1:12

SHEET 1 OF 1

NESTED INSTRUMENT TABLE

are a space saving alternative with a unique design that requires only 8sq.of floor space but yields 31 sq.ft. of usable workspace. They can be order individually or as complete set.



PO#: _____ DATE: _____	 BLICKMAN	NAME	DATE
APPROVAL: _____		DRAWN	AP 04/10/2006
TITLE: _____		ENG APPR.	
		MATERIAL	N/A
		FINISH	N/A
		APPRX WEIGHT = 38.966 lbs	
DO NOT SCALE DRAWING			

BLICKMAN HEALTH INDUSTRIES, INC.
LODI, N.J.

TITLE:

Nested Table
7822SS

SIZE
A

DWG. NO.
0137822000

REV
PF

SCALE: 1:20





SHEET 1 OF 1

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF BLICKMAN HEALTH INDUSTRIES. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF BLICKMAN HEALTH INDUSTRIES IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:
DIMENSIONS ARE IN INCHES
BEND RADIUS = MTL. THK.
TOLERANCES:
FRACTION: $\pm 1/32"$
DECIMAL: $\pm .031"$

NESTED INSTRUMENT TABLE

are a space saving alternative with a unique design that requires only 8sq.of floor space but yields 31 sq.ft. of usable workspace. They can be order individually or as complete set.



PO#: _____ DATE: _____

APPROVAL: _____

TITLE: _____



PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF BLICKMAN HEALTH INDUSTRIES. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF BLICKMAN HEALTH INDUSTRIES IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:
DIMENSIONS ARE IN INCHES
BEND RADIUS = MTL. THK.
TOLERANCES:
FRACTION: $\pm 1/32"$
DECIMAL: $\pm .031"$

	NAME	DATE
DRAWN	AP	04/10/2006
ENG APPR.		
MATERIAL N/A		
FINISH N/A		
APPRX WEIGHT = 38.966 lbs		
DO NOT SCALE DRAWING		

BLICKMAN HEALTH INDUSTRIES, INC.
LODI, N.J.

TITLE:





Nested Table
7824SS

SIZE	DWG. NO.	REV
A	0137824000	PE

SCALE: 1:20 SHEET 1 OF 1

NESTED INSTRUMENT TABLE

are a space saving alternative with a unique design that requires only 8sq.of floor space but yields 31 sq.ft. of usable workspace. They can be order individually or as complete set.



PO#: _____ DATE: _____

APPROVAL: _____

TITLE: _____



PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF BLICKMAN HEALTH INDUSTRIES. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF BLICKMAN HEALTH INDUSTRIES IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:
DIMENSIONS ARE IN INCHES
BEND RADIUS = MTL. THK.
TOLERANCES:
FRACTION: $\pm 1/32"$
DECIMAL: $\pm .031"$

	NAME	DATE
DRAWN	AP	04/10/2006
ENG APPR.		
MATERIAL N/A		
FINISH N/A		
APPRX WEIGHT = 38.966 lbs		
DO NOT SCALE DRAWING		

BLICKMAN HEALTH INDUSTRIES, INC.
LODI, N.J.

TITLE:

Nested Table
7827SS

SIZE A	DWG. NO. 0137827000	REV PE
------------------	------------------------	-----------

SCALE: 1:14

SHEET 1 OF 1