


1.3 mL capacity – small




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134010	Non-sterile Non-sterile	500 pcs/dispenser box, 8 boxes/case 500 pcs/dispenser box	4000 500	<ul style="list-style-type: none"> Capacity: 1.3 mL Drops/mL: 45 Bulb draw: 1 mL Total length: 51 mm (2.0") Stem diameter: 1.15 mm (0.05") Tip diameter: 1.15 mm (0.05")

1.5 mL capacity – extended tip – small bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134020	Non-sterile Non-sterile	400 pcs/dispenser box, 8 boxes/case 400 pcs/dispenser box	3200 400	<ul style="list-style-type: none"> Capacity: 1.5 mL Drops/mL: 50 Bulb draw: 1 mL Total length: 104 mm (4.1") Stem diameter: 4.5 mm (0.18") Tip diameter: 0.9 mm (0.04")

5 mL capacity – extended tip – large bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134070	Non-sterile Non-sterile	250 pcs/dispenser box, 8 boxes/case 250 pcs/dispenser box	2000 250	<ul style="list-style-type: none"> Capacity: 5 mL Drops/mL: 50 Bulb draw: 3.2 mL Total length: 153 mm (6.0") Stem diameter: 4.6 mm (0.18") Tip diameter: 0.85 mm (0.03")

5.8 mL capacity – standard




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134050	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 5.8 mL Drops/mL: 60 Bulb draw: 3.4 mL Total length: 147 mm (5.8") Stem diameter: 4.6 mm (0.18") Tip diameter: 0.85 mm (0.03")
134050-S20	Sterile	20 pcs/pack, 20 packs/case	400	

5.8 mL capacity – extended tip – standard



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134060	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 5.8 mL Drops/mL: 50 Bulb draw: 3.4 mL Total length: 157 mm (6.2") Stem diameter: 4.6 mm (0.18") Tip diameter: 0.85 mm (0.03")

8.7 mL capacity – large bulb



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
134090	Non-sterile Non-sterile	400 pcs/dispenser box, 10 boxes/case 400 pcs/dispenser box	4000 400	<ul style="list-style-type: none"> Capacity: 8.7 mL Drops/mL: 45 Bulb draw: 4.8 mL Total length: 147 mm (5.8") Stem diameter: 6.1 mm (0.24") Tip diameter: 1.4 mm (0.06")