



1.5 mL capacity – graduated to 0.3 mL – pediatric




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
136036	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 1.5 mL Drops/mL: 21 Bulb draw: 0.9 mL Total length: 116 mm (4.6") Stem diameter: 3.2 mm (0.13") Tip diameter: 3.2 mm (0.13") Graduated up to: 0.3 mL Graduation increments: 0.1 mL

3 mL capacity – graduated to 1 mL – small bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
137035	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 3 mL Drops/mL: 25 Bulb draw: 2.2 mL Total length: 140 mm (5.5") Stem diameter: 4.4 mm (0.17") Tip diameter: 2.4 mm (0.09") Graduated up to: 1 mL Graduation increments: 0.25 mL
137135	Sterile Sterile	Individually wrapped, 100/bag, 4 bags/case Individually wrapped, 100/bag	400 100	

5 mL capacity – graduated to 1 mL – large bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
137030	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 5 mL Drops/mL: 24 Bulb draw: 3.4 mL Total length: 145 mm (5.7") Stem diameter: 4.8 mm (0.19") Tip diameter: 2.7 mm (0.11") Graduated up to: 1 mL Graduation increments: 0.25 mL
137038	Sterile Sterile	Individually wrapped, 100/bag, 4 bags/case Individually wrapped, 100/bag	400 100	
137138	Sterile Sterile	10 pcs/pack, 10 packs/bag, 4 bags/case 10 pcs/pack, 10 packs/bag	400 100	
137238	Sterile	20 pcs/pack, 20 packs/case	400	

5 mL capacity – graduated to 1 mL – large bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
137010	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 5 mL Drops/mL: 25 Bulb draw: 3.4 mL Total length: 150 mm (5.9") Stem diameter: 4.8 mm (0.19") Tip diameter: 2.7 mm (0.11") Graduated up to: 1 mL Graduation increments: 0.25 mL
137018 NEW!	Sterile	Individually wrapped peel packs, 400 peel packs/case	400	

5 mL capacity – graduated to 2.25 mL




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
138070	Non-sterile Non-sterile	500 pcs/dispenser box, 10 boxes/case 500 pcs/dispenser box	5000 500	<ul style="list-style-type: none"> Capacity: 5 mL Drops/mL: 21 Bulb draw: 1.8 mL Total length: 150 mm (5.9") Stem diameter: 7.2 mm (0.28") Tip diameter: 3 mm (0.12") Graduated up to: 2.25 mL Graduation increments: 0.25 mL

7 mL capacity – graduated to 3 mL – large bulb



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
135030	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 7 mL Total length: 160 mm (6.3") Drops/mL: 20 Stem diameter: 7.9 mm (0.31") Bulb draw: 3.2 mL Tip diameter: 3 mm (0.12") Graduated up to: 3 mL Graduation increments: 0.5 mL
	Non-sterile	500 pcs/dispenser box	500	
135038	Sterile	Individually wrapped, 100/bag, 4 bags/case	400	
	Sterile	Individually wrapped, 100/bag	100	
135138	Sterile	10 pcs/pack, 10 packs/bag, 4 bags/case	400	
	Sterile	10 pcs/pack, 10 packs/bag	100	
135238	Sterile	20 pcs/pack, 20 packs/case	400	

7.5 mL capacity – graduated to 3 mL – large bulb



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
135010	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 7.5 mL Total length: 148 mm (5.8") Drops/mL: 21 Stem diameter: 7.5 mm (0.30") Bulb draw: 3.5 mL Tip diameter: 3.2 mm (0.13") Graduated up to: 3 mL Graduation increments: 0.5 mL
	Non-sterile	500 pcs/dispenser box	500	